



Branan Medical Corporation

10015 Muirlands Road, Suite C • Irvine, CA 92618 • Tel: (949) 598-7166 • Fax: (949) 598-7167

BMC LAUNCHES ORATECT[®] III ORAL FLUID DRUG SCREEN DEVICE

November 1, 2006

Branan Medical Corporation is pleased to announce the launch of the new and improved Oratect[®] III Oral Fluid Drug Screen Device. Oratect[®] III features the new blue lines in the test windows to indicate when a sufficient amount of saliva is obtained to activate the test. With Oratect[®] III the test results are stable up to thirty minutes and has increased sensitivity.

Oratect[®] III provides a uniquely simple 1-step procedure for the accurate detection of 6 drugs of abuse in oral fluids including marijuana (THC), amphetamine, methamphetamine/ecstasy, cocaine, opiates and either PCP or benzodiazepines. Collection and testing entail swabbing the collection pad inside the mouth until blue lines in the test panels start to flow, indicating a sufficient amount of saliva has been collected. The presence of a red line corresponding to a specific drug on the panel indicates a negative result, while the absence of a red line indicates a presumptive positive result. Presumptive positive results are confirmed by sending the collection pad of the Oratect[®] III device in an included vial to a SAMHSA laboratory.

Oratect[®] III offers the perfect solution for Human Resource Managers and Criminal Justice Professionals to minimize sample adulteration and contamination, and eliminate the need to handle a urine specimen by the administrator.

For more information on Branan Medical Corporation or its drug testing products, please visit www.brananmedical.com.

About Branan Medical Corporation

Branan Medical Corporation, located in Irvine, California, is a manufacturer and marketer of rapid, on-site tests for drugs of abuse. Its products include adulteration tests, CLIA-waived urine drug screens, oral fluid drug and alcohol screens, and drug test controls.