



Branan Medical Corporation

• 10015 Muirlands Road, Suite C • Irvine, CA 92618 • Tel: (949) 598-7166 • Fax: (949) 598-7167

ToxCup[®] CLIA-Waived Press Release

ToxCup[®] Drug Screen Cup is now CLIA-Waived!

Branan Medical Corporation's **ToxCup[®] Drug Screen Cup** has achieved **CLIA-Waived status**. ToxCup[®] Drug Screen Cup is a one-step test for the rapid detection of up to seven drugs of abuse integrating sample collection, drug testing and shipping for confirmation in a single container. "CLIA-Waived" status indicates that the test is easy enough to be used by laboratory and non-laboratory personnel at point of care sites such as treatment centers, occupational medicine departments, hospital laboratories, physician offices, probation/parole agencies, prisons and the workplace.

Branan Medical Corporation is the only company to provide a CLIA-Waived test cup for abused drugs including: Amphetamine, Barbiturates, Benzodiazepines, Cocaine, MDMA, Methamphetamine, Methadone, Opiates, Oxycodone, PCP and THC. ToxCup[®] Drug Screen Cup is easy to use in any facility. No instrumentation is required. Twist the test lid onto the specimen cup after collection. Then lean the cup on its side to activate testing. The test is completed in 5 minutes or less and includes an internal procedural control.

ToxCup[®], supplied 25 devices per box, is available in various 5-drug, 6-drug, 7-drug configurations depending on the needs of your facility. ToxCup[®] has an 18-month shelf life and can be stored at room temperature (15°-30°C).

For more information, please call Branan Medical Corporation at 949-598-7166 or e-mail us at info@brananmedical.com.