



## Branran Medical Corporation

---

• 10015 Muirlands Road, Suite C • Irvine, CA 92618 • Tel: (949) 598-7166 • Fax: (949) 598-7167

### **Fastect<sup>®</sup> II Drug Screen Dipstick Test is now CLIA-Waived!**

Branran Medical Corporation's Fastect<sup>®</sup> II Drug Screen Dipstick Test, an on-site device for the simultaneous detection of 1 to 4 commonly abused prescription or illegal drugs in urine, has achieved **CLIA-Waived status**. "CLIA-Waived" status indicates that the test is easy enough to be used by laboratory and non-laboratory point of care personnel in treatment centers, occupational medicine departments, hospital labs, physician offices, probation/parole agencies, prisons and workplace test sites.

Branran Medical Corporation is the only company to provide a CLIA-Waived test device for the most commonly abused drugs including: Marijuana, Cocaine, Opiates, Methamphetamine, Amphetamine, Benzodiazepines, Barbiturates, Oxycodone, MDMA, PCP, Methadone, and Tricyclic Antidepressants.

Fastect<sup>®</sup> II Drug Screen Dipstick Test is easy to use in any facility. No instrumentation is required. The device is dipped in the sample for 10 seconds and removed. The test is completed in 5 minutes or less and includes internal procedural controls.

Fastect<sup>®</sup> II, supplied 50 devices per box, is available in single or multiple configurations depending on the needs of your facility. It has an 18-month shelf life and can be stored at room temperature.

For more information, please call Branran Medical Corporation at 949-598-7166 or e-mail us at [info@brananmedical.com](mailto:info@brananmedical.com).