



Branran Medical Corporation

• 10015 Muirlands Road, Suite C • Irvine, CA 92618 • Tel: (949) 598-7166 • Fax: (949) 598-7167

ORATECT® II HM12 Press Release

Irvine, CA – April 3, 2006 - Branran Medical Corporation announces the introduction of **Oratect® II HM12**, the newest addition to the Oratect II product line. HM12 simultaneously tests for six drugs in oral fluid: marijuana (THC), cocaine, opiates, amphetamine, methamphetamine (including MDMA) and benzodiazepines.

Oratect II is a unique oral fluid drug screen test that incorporates the collection and testing processes into a single device. It allows a simplified one-step procedure while decreasing the time required to obtain a test result to 7 minutes from start to finish.

The new Oratect II product line features:

- Two device configurations: HM11 with PCP; HM12 with benzodiazepines.
- Indication to show sufficient amount of saliva is collected.
- Increased drug test sensitivities (with parent THC at 40ng/ml) to provide a longer window of detection.
- Vial included for further testing allowing confirmation of the onsite result without drawing of another sample.

This product is supplied with 25 devices per box and is available through your local distributor. It is intended for forensic use.

Branran Medical Corp is a leading U.S. manufacturer of drugs of abuse testing kits located in Irvine, CA. For more information, call toll free: 866-468-3287, worldwide: 949-598-7166 or send an email message to info@brananmedical.com.